

ED Agency Urine Pregnancy Test

Name: _____ Date: _____ Agency: _____

1. Which specimen contains the highest concentration of HCG and therefore is the best sample?
 - a) The first morning urine sample
 - b) A mid-day specimen
 - c) An evening specimen
 - d) Any time of day specimen
2. What is the minimum amount of time you should wait after applying the sample before deciding that a result is negative?
 - a) 1 minute
 - b) 2 minutes
 - c) 3 minutes
 - d) 4 minutes
3. You should read the results in how many minutes?
 - a) 3 minutes
 - b) 15 minutes
 - c) 30 minutes
 - d) 60 minutes
4. After interpretation of the test, the results should be documented in which of the following places:
 - a) In the computer
 - b) On the manual log with a patient sticker
 - c) All of the above
5. Which of the following represents an invalid result?
 - a) A line in the "C" window, no line in the "T" window
 - b) No line in the "C" window, no line in the "T" window
 - c) A line in both windows